



MEDICO

TECHNOLOGIES



Dental Instruments

Catalog **2021**



EXTRACTION FORCEPS

APPLICATION CHART FOR EXTRACTION FORCEPS





Fig	Application
1	Upper centrals, canines and roots
2	Upper laterals
7	Upper premolars and roots
13	Lower premolars
17	Upper molars (Right)
18	Upper molars (Left)
19	Upper wisdoms
22	Lower molars (Hawk's bill)
29	Upper roots
29S	Upper roots
30	Upper roots
37	Children's upper incisors and canines
39	Children's upper molars
44N	Upper roots
51	Upper wisdoms (Bayonet)
67	Upper wisdoms
73	Lower molars (Hawk's bill)
73S	Lower molars (Hawk's bill)
74	Lower roots
74N	Small lower roots and anteriors
75	Lower premolars
76	Upper roots (Read)
76S	Upper roots (Read)
76N	Small upper roots (Read)

Fig	Application
79	Lower wisdoms
86	Lower molars (Cowhorn)
87	Lower molars (Cowhorn)
89	Upper molars (Right)(Upper Cowhorns)
90	Upper molars (Left)(Upper Cowhorns)
94	Upper molars (Right)
95	Upper molars (Left)
101	Upper molars and roots
107	Upper canines
113	Upper roots (Martin)
123	Children's lower anteriors, premolars, roots
136	Upper premolars and roots (Guy's)
137	Upper anteriors and roots (Guy's)
138	Children's upper anteriors, premolars, roots
147	Upper roots (Gibb)
157	Children's upper molars
158	Children's upper molars
159	Children's upper premolars
160	Children's lower molars (Hawk's bill)
161	Children's lower molars (Hawk's bill)
162	Children's lower anteriors and roots
163	Children's upper anteriors
3	Rongeurs forceps



HANDLE OPTIONS

SMALL MODEL



LARGE MODEL



FINGER-SHAPED



DIAMOND TIPPED





EXTRACTION FORCEPS



MT-450-01 Fig. 01
Upper centrals and canines, wide



MT-450-02 Fig. 02
Upper laterals and canines



MT-450-03 Fig. 03
Upper crowded incisors & canines



MT-450-04 Fig. 04
Lower incisors and canines



MT-450-07 Fig. 07
Upper Premolars



MT-450-08 Fig. 08
Lower premolars either side



MT-450-13 Fig. 13
Lower premolars either side



MT-450-13S Fig. 13S
Lower premolars, for children



MT-450-17 Fig. 17
Upper molars, right




EXTRACTION FORCEPS



MT-450-18 **Fig. 18**
Upper molars left



MT-450-18A **Fig. 18A**
Upper molars either side



MT-450-18B **Fig. 18B**
Upper molars



MT-450-19 **Fig. 19**
Upper third molars & upper wisdoms either side



MT-450-20 **Fig. 20**
Lower third molars



MT-450-21 **Fig. 21**
Lower molars



MT-450-22 **Fig. 22**
Lower molars



MT-450-22S **Fig. 22S**
Lower molars, for children



MT-450-23 **Fig. 23**
Lower molars, right



EXTRACTION FORCEPS



MT-450-24 Fig. 24
Lower molars, left



MT-450-29 Fig. 29
Upper roots and incisors



MT-450-29N Fig. 29N
Root of Upper Incisors



MT-450-29S Fig. 29S
Upper roots, for children



MT-450-30 Fig. 30
Upper roots



MT-450-31 Fig. 31
Lower roots & incisors with
Serrated tips



MT-450-32 Fig. 32
Lower molars



MT-450-33 Fig. 33
Lower roots



MT-450-33A Fig. 33A
Lower roots



EXTRACTION FORCEPS



MT-450-33C Fig. 33C
Lower roots




MT-450-33L Fig. 33L
Lower roots



MT-450-33M Fig. 33M
Lower roots with long beak



MT-450-34 Fig. 34
Incisors and Bicuspid,
deep gripping



MT-450-34N Fig. 34N
Incisors and Bicuspid,
deep gripping



MT-450-35 Fig. 35
Incisors and Bicuspid,
deep gripping



MT-450-35N Fig. 35N
Incisors and Bicuspid,
deep gripping



MT-450-36 Fig. 36
Incisors and Bicuspid,
deep gripping



MT-450-37 Fig. 37
Upper incisors, for children




EXTRACTION FORCEPS




MT-450-38 **Fig. 38**
Lower incisors and
canines, for children




MT-450-39 **Fig. 39**
Upper molars, for children



MT-450-39A **Fig. 39A**
Upper molars, for children



MT-450-39L **Fig. 39L**
Upper molars left, for children



MT-450-39R **Fig. 39R**
Upper molars right, for children



MT-450-40 **Fig. 40**
Lower molars, for children



MT-450-41 **Fig. 41**
Upper roots



MT-450-44 **Fig. 44**
Upper roots



MT-450-44F **Fig. 44F**
Upper roots



EXTRACTION FORCEPS




MT-450-45 Fig. 45
Lower roots



MT-450-45F Fig. 45F
Lower roots



MT-450-46 Fig. 46
Lower roots




MT-450-46L Fig. 46L
Lower roots



MT-450-46LX Fig. 46LX
for very fine roots



MT-450-46N Fig. 46N
Bicuspid, deep gripping



MT-450-49 Fig. 49
Upper roots




MT-450-50 Fig. 50
Upper roots



MT-450-51 Fig. 51
Upper roots



EXTRACTION FORCEPS




MT-450-51A Fig. 51A
Upper roots



MT-450-51C Fig. 51C
Upper roots



MT-450-51L Fig. 51L
Upper roots



MT-450-51M Fig. 51M
Especially long beaks



MT-450-52 Fig. 52
Upper roots



MT-450-52A Fig. 52A
Upper teeth & Resection



MT-450-53 Fig. 53
Upper roots



MT-450-54 Fig. 54
For separating upper molars




MT-450-55 Fig. 55
For separating upper molars




EXTRACTION FORCEPS



MT-450-56 **Fig. 56**
For separating lower molars



MT-450-65L **Fig. 65L**
Upper molars, left




MT-450-65R **Fig. 65R**
Upper molars, right




MT-450-67 **Fig. 67**
Upper third molars and
upper wisdoms either side



MT-450-67A **Fig. 67A**
Upper third molars and upper
wisdoms either side



MT-450-67LX **Fig. 67LX**
Upper third molar, left



MT-450-67RX **Fig. 67RX**
Upper third molar, right



MT-450-67N **Fig. 67N**
Third molar, deep-gripping



MT-450-68 **Fig. 68**
Lower roots



EXTRACTION FORCEPS



MT-450-73 Fig. 73
Lower molars



MT-450-73L Fig. 73L
Lower molars



MT-450-73S Fig. 73S
Lower molars



MT-450-74 Fig. 74
Lower roots and incisors



MT-450-74D Fig. 74D
Lower roots



MT-450-74M Fig. 74M
Lower roots and incisors



MT-450-74N Fig. 74N
Lower roots and incisors



MT-450-75 Fig. 75
Lower incisors and premolars




MT-450-75A Fig. 75A
Lower incisors and premolars




EXTRACTION FORCEPS




MT-450-76 **Fig. 76**
Upper roots



MT-450-76N **Fig. 76N**
Upper roots




MT-450-76S **Fig. 76S**
Upper roots



MT-450-77 **Fig. 77**
Upper roots



MT-450-79 **Fig. 79**
Lower third molars and
lower wisdoms either side



MT-450-79A **Fig. 79A**
Lower third molars



MT-450-79N **Fig. 79N**
Lower Wisdoms



MT-450-86 **Fig. 86**
Lower molars decayed
or broken down crowns



MT-450-86B **Fig. 86B**
Lower molars



EXTRACTION FORCEPS



MT-450-86C **Fig. 86C**
Lower molars decayed or
broken down crowns



MT-450-87 **Fig. 87**
Lower molars decayed or
broken down crowns




MT-450-89 **Fig. 89**
Upper molars, right




MT-450-90 **Fig. 90**
Upper molars, left



MT-450-91 **Fig. 91**
For separating lower molars



MT-450-92 **Fig. 92**
Upper molars, right



MT-450-93 **Fig. 93**
Upper molars, left



MT-450-94 **Fig. 94**
Upper molars, right



MT-450-95 **Fig. 95**
Upper molars, left




EXTRACTION FORCEPS



MT-450-96 **Fig. 96**
Upper molars, right




MT-450-97 **Fig. 97**
Upper molars, left




MT-450-99 1/4 **Fig. 99 1/4**
Lower molars



MT-450-99 **Fig. 99**
Lower molars



MT-450-101 **Fig. 101**
Upper roots




MT-450-101A **Fig. 101A**
Upper roots



MT-450-107 **Fig. 107**
Upper Canines




MT-450-108 **Fig. 108**
Upper roots




MT-450-110 **Fig. 110**
Upper roots



EXTRACTION FORCEPS



MT-450-111 Fig. 111
Upper roots



MT-450-112 Fig. 112
Upper roots



MT-450-113 Fig. 113
Upper roots



MT-450-114 Fig. 114
Lower roots & centrals



MT-450-121 Fig. 121
Upper wisdoms



MT-450-122 Fig. 122
Lower centrals & roots



MT-450-123 Fig. 123
Lower teeth & roots,
for children



MT-450-130 Fig. 130
Upper wisdoms 190 mm




MT-450-133 Fig. 133
for very fine roots



EXTRACTION FORCEPS



MT-450-136 **Fig. 136**
Upper Centrals,
Premolars & roots



MT-450-136H **Fig. 136H**
Upper roots



MT-450-137 **Fig. 137**
Lower centrals & roots



MT-450-138 **Fig. 138**
Upper teeth & roots,
for children




MT-450-139 **Fig. 139**
Upper roots




MT-450-143 **Fig. 143**
Lower roots & centrals



MT-450-144 **Fig. 144**
Lower roots & centrals



MT-450-145 **Fig. 145**
Lower roots



MT-450-146 **Fig. 146**
Root fragments,
deep gripping




EXTRACTION FORCEPS




MT-450-147 Fig. 147
Upper roots




MT-450-148 Fig. 148
Lower molars



MT-450-151 Fig. 151
for very fine roots



MT-450-153 Fig. 153
Root fragments,
deep gripping



MT-450-155 Fig. 155
Lower roots



MT-450-157 Fig. 157
Upper molars,
for children



MT-450-158 Fig. 158
Upper premolars & molars,
for children



MT-450-159 Fig. 159
Upper premolars,
for children



MT-450-160 Fig. 160
Lower molars,
for children



EXTRACTION FORCEPS



MT-450-161 **Fig. 161**
lower molars, for children



MT-450-162 **Fig. 162**
Lower roots & centrals,
for children



MT-450-163 **Fig. 163**
Upper centrals, for children



MT-450-166 **Fig. 166**
Upper universal application



MT-450-167 **Fig. 167**
Lower teeth, universal pattern



MT-450-168 **Fig. 168**
Upper central & lateral teeth



MT-450-169 **Fig. 169**
Premolars, lower canines
& all roots




MT-450-171 **Fig. 171**
Wisdoms & molar roots



MT-450-173 **Fig. 173**
Lingually inclined roots and
wisdom teeth



EXTRACTION FORCEPS



MT-450-180 **Fig. 180**
Special for upper centrals and
canines deep gripping




MT-450-181 **Fig. 181**
Special for upper premolars,
deep gripping



MT-450-182 **Fig. 182**
Lower centrals canines and
premolars deep gripping



MT-450-185 **Fig. 185**
Lower roots



MT-450-187 **Fig. 187**
Lower roots fragments



MT-450-189 **Fig. 189**
Fragments upper roots



MT-450-190 **Fig. 190**
Fragments upper roots



MT-450-221 **Fig. 221**
Universal upper roots 150 mm



MT-450-223 **Fig. 223**
Universal for lowr roots



EXTRACTION FORCEPS



MT-450-350 **Fig. 350**
Lower molars & wisdoms left



MT-450-351 **Fig. 351**
Lower molars & wisdoms right

MEDICO TECHNOLOGIES



wazirabad road sialkot- 51310 pakistan



0092 523 571 255



0092 523 571 255



info@medico-medtech.com



<http://www.medico-medtech.com>