



MEDICO

TECHNOLOGIES



Electrosurgical Instruments

Catalog **2021**



MONOPOLAR CABLES



Silicone Monopolar Cables 10" (3.0m)

MT - RB 19.230

Instrument Connector: Screw Type
Generator Connector: Standard Bovie
8.0mm Plug Compatibility: Valleylab,
Conmed, Bovie



Silicone Monopolar Endoscopic Cable 10" (3.0m)

MT - RB 19.240

Instrument Connector: 4.0mm Female
Plug Generator Connector: Standard
Bovie 8.0mm Plug Compatibility:
Valleylab, Conmed, Bovie



Silicone Monopolar Endoscopic Cable 10" (3.0m)

MT - RB 19.250

Instrument Connector: 4.0mm Female
Plug Generator Connector: Standard
4.0mm Banana Plug Compatibility:
Valleylab, Conmed & All ES Generators



Silicone Monopolar Cables 10" (3.0m)

MT - RB 19.260

Instrument Connector: 4.0mm Male
Fitting Generator Connector: Standard
4.0mm Banana Plug Compatibility:
Valleylab, Conmed & All ES Generators



Silicone Diathermy Cable 10" (3.0m)

MT - RB 19.270

Instrument Connector: 4.8mm Male
Fitting Generator Connector: Standard
Bovie 8.0mm Plug Compatibility:
Valleylab, Conmed, Bovie



Silicone Diathermy Cable 10" (3.0m)

MT - RB 19.260

Instrument Connector: 4.8mm Male
Fitting Generator Connector: Standard
4.0mm Banana Plug Compatibility:
Valleylab, Conmed & All ES Generators

MEDICO TECHNOLOGIES



wazirabad road sialkot- 51310 pakistan



0092 523 571 255



0092 523 571 255



info@medico-medtech.com



<http://www.medico-medtech.com>